

B. Huquet¹, F. Barillet³, M. Brochard¹,

G. Miranda³, H. Larroque³, O. Leray⁴,

PhénoFinlait Cheap and large scale analysis methods to quantify milk fatty acids and proteins content

Context

- Evolution of milk quality concept: nutritive > sanitary clean > healthy fine composition
- Lack of analysis methods for milk fatty acid (FA) and individual protein (IP) usable at high scale


- Implication of every partners of french dairy industry, 1,500 commercial
- First objectives: develop analysis method to determine milk FA profile, create reference analysis method to determine milk IP profile

The french answer: PhénoFinlait research program

farms, 3 species (cow, ewe, goat) and 7 breeds

General principle of a cheap and large scale analysis method

Reference analysis



F. Faucon-Lahalle^{1,2},

I. Palhière³, J.L. Peyraud³,

P. Martin³ and M. Ferrand¹

¹Institut de l'Elevage;

²CNIEL; ³INRA;

⁴Actilait

MIR spectra reflects biochemical composition as absorbance reflects biochemical bonds of molecules

MIR* spectra Routinely obtained by milk recording laboratory for fat and protein percentage measurements/

Statistical Model 3

Calibration **Equations**

PhénoFinlait program, collection of 20,000 MIR spectra/week/lab

MIR spectra calibration

with control milks

Estimation of fine milk composition (FA, IP)




* MIR = Mid Infra-Red

For fatty acids

- Gas chromatography: 400 analysis for cows, 200 for ewes and 200 for goats
- a) Selection of wavelengths with genetic algorithms (Leardi R. 1998, Ferrand 2009) b) Partial Least Square (PLS) regression


(coefficients are applied to wavelengths)

Estimations usable for 15 to 20 FA



For individual proteins

Establishment of the reference method:



- Calculation with 67 cow milk samples Partial Least Square (PLS) regression => coefficients are applied to wavelengths
- Too many variants (genetic, postranscriptional...) to have good estimations of IP Necessity to analyze more milk samples (14,000 in the PhenoFinlait program)

